

SPECIFICATIONS

MicroActive Share™

INSTRUMENT CONNECTION TO MICROACTIVE SHARE CLOUD

Micromeritics 3Flex or ASAP 2020 MA Share instrument Data Controller/Module with active Internet connection for upload of data to MicroActive Share Cloud.

FOR USER CONNECTION TO MICROACTIVE SHARE CLOUD

- Internet connected device with installed browser (SmartPhone, Tablet, Desktop/laptop)
- Users must have a valid and active email address

BROWSER COMPATIBILITY

Mac OS X	<ul style="list-style-type: none">• Safari 7 .x• Google Chrome• Firefox 4 or later
Windows 7	<ul style="list-style-type: none">• Internet Explorer 10• Internet Explorer 11• Google Chrome• Firefox 4 or later
Windows 8	<ul style="list-style-type: none">• Internet Explorer 10• Internet Explorer 11• Google Chrome• Firefox 4 or later
Windows 10	<ul style="list-style-type: none">• Internet Explorer 10• Internet Explorer 11• Google Chrome• Firefox 4 or later

AVAILABLE REPORTS

Available Reports	<ul style="list-style-type: none">• Isotherm,• BET Surface Area,• t-Plot• BJH Adsorption and Desorption• Horvath-Kawazoe,• NLDFT
-------------------	---

SPECIFICATIONS

MicroActive Share

SECURITY PROTOCOLS

Security Protocols

- ISO 27002, ISO 27001
- SOC 1, SOC2, SOC 3
- Safe Harbor, SSAE16, PCI-DSS and CPS

SITE STRUCTURE

Site Structure

A site is created when an MA Share Instrument Owner is created. Administrator can create a collection of members that have an assigned level of access rights.

USER HIERARCHY

Site Owner

Owner/purchaser of MA Share instrument- Full rights

Administrator

Ability to add or delete members from site, capable of promoting and demoting rights of members within site.

Can create public or private groups

Member

Access to data at assigned level (read or read/write)

Can create public or private groups

RESULT SETS

Result Sets

- **My Samples:** Data stored under assigned user names
- **Site Samples:** Data stored under a site
- **Public Samples:** Samples that have been public to site members

MICROACTIVE SHARE COMMUNITY

Public discussion forum for information exchange. User can post to public forum and receive responses from the members within the MicroActive community. Members must have a MA Share account to contribute to MA Share Community. Post can directly be loaded to social media.

**Due to continuous improvements, specifications are subject to change without notice.*